

# 結案報告 系統功能展示

---

日期：2020/06/03

單位：工程行銷－小組（外包 陳淳濤）

# 大綱

---

- 一. 帳號角色介紹
- 二. 功能概要
- 三. 操作流程

# 帳號角色介紹

---

一. 病人

二. 醫生

# 有哪些使用者

---

## 1. 病人

- 一般網路使用者
- 透過登入頁面註冊獲取帳號

## 2. 醫生

- 系統內部使用者
- 透過管理員給予新的帳號

# 功能概要

---

- 本專案開發了那些新功能

# 本專案開發了那些新功能

---

## 1. 系統登入

- 填寫病歷資料

## 2. 病人

- 醫師留言通知
- 瀏覽個人病歷資料
- 填寫個人量測記錄
- 瀏覽個人歷史量測記錄

## 3. 醫生

- 查詢病人病歷資料
- 查詢病人歷史量測記錄與給予評估意見
- 設定系統量測記錄警訊值



# 操作流程

---

- 一. 病人部分
- 二. 醫生部分

# 病人部分

---

## Part 1：獲取帳號

1. 註冊帳號
2. 填寫病歷資料

## Part 2：病人功能操作

1. 量測資料填寫
2. 檢視歷史量測資料
3. 查看醫師建言



# PART I : 獲取帳號 註冊帳號

**Health Care Login & Registration Form**

THIS IS A PROTOTYPE and registration form for HEALTH CARE CARE

**Login to our Website**

ENTER LOGIN ID AND PASSWORD TO LOG IN

ACCOUNT

LOGIN

PASSWORD

REGISTER

[I forgot my password](#)

[or login with](#)

[+ add social](#)

**Sign up now!**

It is the right way to get patient access

Name

Birth Date

Birth Date


Gender

Phone Number

Address

[I've had an account](#)

[Sign Up](#)



**Sign up now**

Fill in the form below to get instant access

Name

Birth Date

Email

Password

Confirm Password

Phone Number

Gender

Address

[I've had an account](#)

[Sign Up](#)

# PART I : 獲取帳號 填寫病歷資料



# PART I : 獲取帳號 填寫病歷資料



Initial login page with a large red arrow pointing to the next step.



Initial Information (初始問診) form. It includes fields for patient name, date of birth, gender, and a list of symptoms.



Detailed History (二、詳細病史) form. It includes a table for symptoms and a section for medical history.

症狀	發	熱	寒	汗
發熱	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
畏寒	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
自汗	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
盜汗	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
手足心熱	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
口乾	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



Personal Information (三、個人基本資料) form. It includes fields for name, date of birth, gender, and a list of symptoms.

# PART I : 獲取帳號 完成註冊

The screenshot shows a web browser window with a dark header bar containing navigation links in Chinese. The main content area has a blue background with the title "Your Information". A white form box is centered, containing the following fields:

- Name: 潘志忠
- Phone: 0933066279
- Birth date: 2023-04-08
- Age: 8
- Gender: M
- Address: TW Taichung

The bottom of the page features a dark footer with contact information and a company description.

**Work Contact Info**


- Address: 14 Rue de l'Indépendance, Lahore
- Phone: 8122-7761894
- Email: support@rethink.com

**About the company**


Our mission is to provide a secure and reliable platform for businesses to manage their data and operations. We are committed to providing the best service and support to our customers.



## PART 2 : 病人功能操作 功能導覽



## PART 2 : 病人功能操作 量測資料填寫(支援兩種填寫模式)



### MedicX 4 Health Care 填寫個人量測記錄表 單

這是一張“量測記錄表單” for Health Care Clinic.

表格模式

表單模式



## PART 2 : 病人功能操作 量測資料填寫(表格模式)

備註：  
BMI值系統會自動計算做帶入

The screenshot shows the MedoX 4 Health Care interface. At the top, it says 'MedoX 4 Health Care 病人個人量測記錄表'. Below this, there's a section for '量測資料填寫' (Measurement Data Entry). The form includes fields for '項目' (Item), '量測數值' (Measurement Value), and '量測時間' (Measurement Time). The items listed are: 身高 (Height), 體重 (Weight), BMI (BMI), 體溫 (Temperature), 脈搏/心跳 (Pulse/Heart Rate), 血氧 (SpO2), 血糖 (Blood Sugar), and 血壓 (Blood Pressure). Each item has a corresponding input field for the value and a date field for the time. The BMI field is highlighted with a red box, indicating it is calculated automatically.

量測記錄表格

請填寫以下量測資料

項目	量測數值	量測時間
身高	請填寫身高(cm)	2020/6/2
體重	請填寫體重(kg)	2020/6/2
BMI BMI參考公式: 體重(公斤) / 身高平方(公尺平方)		2020/6/2
體溫 (標準值: 25 ~ 40)	請填寫體溫	2020/6/2
脈搏/心跳 (標準值: 100 ~ 200)	請填寫脈搏或心跳	2020/6/2
血氧 (標準值: 500 ~ 1000)	請填寫血氧	2020/6/2
血糖 (標準值: 20 ~ 80)	請填寫血糖	2020/6/2
血壓(收縮壓/舒張壓) 收縮 (標準值: 50 ~ 100) 舒張 (標準值: 120 ~ 200)	請填寫血壓	2020/6/2

送出量測資料

## PART 2 : 病人功能操作 量測資料填寫(表單模式)

備註：  
BMI值系統會自動計算做帶入

The screenshot shows a web application interface for 'Medick & Health Care'. A modal form titled '量測資料' (Measurement Data) is displayed in the center. The form contains several input fields for patient information and vital signs, including date, height, weight, temperature, pulse, blood oxygen, blood sugar, and blood pressure. The background is a blue gradient with a faint image of a person.

**量測記錄表單**

請填寫以下量測資料

**量測日期**  
2020/6/2

**身高**  
請填寫身高(cm)

**體重**  
請填寫體重(kg)

**BMI**  
BMI參考公式: 體重(公斤) / 身高平方(公尺平方)

**體溫 (標準值: 25 ~ 40)**  
請填寫體溫

**脈搏/心跳 (標準值: 100 ~ 200)**  
請填寫脈搏或心跳

**血氧 (標準值: 500 ~ 1000)**  
請填寫血氧


**血糖 (標準值: 20 ~ 80)**  
請填寫血糖

**血壓(收縮壓/舒張壓) 收縮 (標準值: 50 ~ 100) 舒張 (標準值: 120 ~ 200)**  
請填寫血壓

送出量測資料

## PART 2 : 病人功能操作 檢視個人歷史量測資料

備註：  
量測時間不會顯示在表中  
請點選【顯示量測時間】觀看



我的歷史量測資料

182.235.107.91/Patient/PatientMeasurementDataView.aspx

	No.	量測表單編號	身高	體重	BMI	體溫	心跳脈搏	血氧	血糖	血壓
顯示量測時間	1	6	174	80	26.424	36.5	100	50.78	68.735	95/130
顯示量測時間	2	7	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	3	8	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	4	9	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	5	10	174	90	29.727	36.5	90	100.78	90.735	75/100
顯示量測時間	6	11	180	75	23.148	40	100	50.78	90.735	95/130

表單填寫時間: 2020/3/5

Health Care © 2017

Milaad St Block B  
Faisal Town, Lahore  
0322 - 7561001  
support@medicx4.com

### About the company

Our mission is to provide quality health services to all in a comfortable and cozy environment. The clinic primarily offers consultations and is not equipped to take emergencies.



## PART 2 : 病人功能操作

### 查看醫師建言

The screenshot displays a web application interface for a patient portal. At the top, a dark navigation bar contains links: Home, 醫師留言通知 (highlighted with a green box), Current Appointment, Take Appointment, Treatment History, 我的病歷資料, 填寫量測資料, 我的歷史量測資料, and Log Out. Below this, a browser window shows the URL 182.235.107.91/Patient/PatientNotifications.aspx. The page title is 醫師留言通知. A message states 'There are no new notifications :)'.

A red box highlights a notification entry:

- 【醫師叮嚀與提醒】
- 發言醫師: AMI
- 留言時間: 2020/6/2 上午 03:33:54
- 留言內容:
- 請小心紅字傷

The footer contains contact information for Health Care © 2017, including the address (Milaad St Block B, Faisal Town, Lahore), phone number (0322 - 7561001), email (support@medicx4.com), and a Facebook link.



# 醫生部分

---

## Part 1：獲取帳號

1. 管理員建立帳號
2. 醫生使用帳號登入

## Part 2：醫生功能操作

1. 查詢與觀看病人病歷資料
2. 查詢觀看病人歷史量測記錄與給予建言
3. 設定系統量測記錄警訊值

# PART I : 獲取帳號 管理員建立帳號

The screenshot displays the Medix4 Health Care system interface. A red box highlights the 'Add Doctor' button in the top navigation bar, which is also highlighted in green. The interface includes a 'Clinic Stats' section on the left, a 'DEPARTMENT INFORMATION' table, and a 'Sign up now' form on the right.

**Clinic Stats**

Total Number of Registered Doctors: 10  
Total Registered Patients: 10  
Total Doctors: 10

**Current Appointments**

**DEPARTMENT INFORMATION**

Department	Description	Number of Doctors
Neurology	Neurology is the study of the nervous system. It includes the brain, spinal cord, and peripheral nerves. Neurologists diagnose and treat disorders of the nervous system.	10
Cardiology	Cardiology is the study of the heart and blood vessels. It includes the diagnosis and treatment of heart disease, high blood pressure, and other conditions.	10
Orthopedics	Orthopedics is the study of the musculoskeletal system. It includes the diagnosis and treatment of injuries and conditions of the bones, joints, and muscles.	10
Internal Medicine	Internal Medicine is the study of the organs and systems of the body. It includes the diagnosis and treatment of a wide range of conditions.	10
Pediatrics	Pediatrics is the study of children and adolescents. It includes the diagnosis and treatment of conditions that affect children.	10
Gynecology	Gynecology is the study of the female reproductive system. It includes the diagnosis and treatment of conditions that affect women.	10

**Sign up now**

ChenChunYun  
06/11/1988  
ccy@gmail.com  
\*\*\*  
\*\*\*  
0900666879  
33000 3000 1  
Neurology PHD IN EVERY FRI ENJOY  
TW,Keelung  
Male Female  
Add



# PART I : 獲取帳號 醫生使用帳號登入

Health Care Login & Registration Form

THIS IS A PROTOTYPE and registration form for Health Care Clinic

Login to our Website

ENTER EMAIL AND PASSWORD TO LOG IN

Account

Email

Password

Remember Me

Login

or login with

Facebook

Sign up now!

It is the right way to get patient access

Name

First Name

Last Name

Phone

Mobile Number


Address

Gender

Male / Female

Register

If you have any queries please feel free to contact us at 011-11111111



Login to our Website

Enter Email and Password to log in:

Account


ccy@gmail.com

Password


...

Login

# PART I : 獲取帳號 登入成功



# PART 2 : 醫生功能操作 功能導覽



## PART 2 : 醫生功能操作 查詢與觀看病人病歷資料


[Home](#) [Patients History](#) [Pending Appointments](#) [Today's Appointments](#) [病人基本病歷資料查詢](#) [設定病人量測資料評估數值](#) [病人量測資料查詢](#) [Log Out](#)

姓名


查詢

	病人姓名	聯繫電話
顯示病歷資料	ABC	61536516
顯示病歷資料	DEF	61536516
顯示病歷資料	XYZ	61536516
顯示病歷資料	neo	0900654888
顯示病歷資料	Test0427	0202000000
顯示病歷資料	陳澤遠	0900666879


Health Care © 2017



Milaad St Block B  
Faisal Town, Lahore




0322 - 7561001



support@medicx4.com

About the company

Our mission is to provide quality health services to all in a comfortable and cozy environment. The clinic primarily offers consultations and is not equipped to take emergencies.





## PART 2 : 醫生功能操作 用姓名查詢

The screenshot displays a web application for medical management. The navigation bar at the top includes links for Home, Patients History, Pending Appointments, Today's Appointments, Patient Basic Medical History Search (highlighted with a green box), Set Patient Measurement Evaluation Value, Patient Measurement Search, and Log Out.

The main content area shows a search form with a text input field labeled '姓名' (Name) and a '查詢' (Search) button. A red box highlights the search input field, and a red arrow points to the search button.

Below the search form, a table displays the search results. The table has columns for '顯示病歷資料' (Display Medical History), '病人姓名' (Patient Name), and '聯絡電話' (Contact Number). The results are as follows:

顯示病歷資料	病人姓名	聯絡電話
顯示病歷資料	ABC	61536516
顯示病歷資料	DEF	61536516
顯示病歷資料	XYZ	61536516
顯示病歷資料	neo	0900654888
顯示病歷資料	Test0427	0202000000
顯示病歷資料	陳淳濤	0900666879

The bottom section of the page contains contact information for the clinic, including the address (Milaad St Block B, Faisal Town, Lahore), phone number (0322 - 7561001), email (support@medicx4.com), and a Facebook link. The footer also includes the text 'Health Care © 2017' and 'About the company'.

## PART 2 : 醫生功能操作 檢視病人病歷資料

備註：  
病人病歷資料不會顯示在表中  
請點選【顯示病歷資料】觀看

姓名:  查詢

顯示病歷資料

	病人姓名	聯絡電話
顯示病歷資料	ABC	61536516
顯示病歷資料	DEF	61536516
顯示病歷資料	XYZ	61536516
顯示病歷資料	neo	0900654888
顯示病歷資料	Test0427	0202000000
顯示病歷資料	陳淳潔	0900666879

初診資料

一、就醫與工作

工務詳

二、西醫

所屬醫院：本國

是否為中國旅行醫：否

中國旅行醫編號

外型記錄

三、家庭病史

有無家族病史：是

家族病史詳細

四、病史

1. 個人內科病史

2. 個人手術史

是否曾患糖尿病：否

家族訂製

0001:00.01

五、個人體格檢查

有無體格檢查：是

體格檢查記錄

六、最近三個月的體格

七、目前之檢驗與藥物

初診資料

一、就醫與工作

工務詳

二、西醫

所屬醫院：本國

是否為中國旅行醫：否

中國旅行醫編號

外型記錄

三、家庭病史

有無家族病史：是

家族病史詳細

四、病史

1. 個人內科病史

2. 個人手術史

是否曾患糖尿病：否

家族訂製

0001:00.01

五、個人體格檢查

有無體格檢查：是

體格檢查記錄

六、最近三個月的體格

七、目前之檢驗與藥物



## PART 2 : 醫生功能操作

### 查詢觀看病人歷史量測記錄與給予建言

[Home](#) [Patients History](#) [Pending Appointments](#) [Today's Appointments](#) [病人基本病歷資料查詢](#) [設定病人量測資料評估數值](#) [病人量測資料查詢](#) [Log Out](#)

姓名  [查詢](#)


顯示資料


	病人姓名	聯絡電話
顯示病人量測資料	ABC	61536516
顯示病人量測資料	DEF	61536516
顯示病人量測資料	XYZ	61536516
顯示病人量測資料	neo	0900654888
顯示病人量測資料	Test0427	0202000000
顯示病人量測資料	陳厚遠	0900666879


**【醫師量測資料意見評估】**  
病人姓名:  
訪視量測資料評估意見

[請按訊息](#) [重新填寫](#)

Healthy Care © 2017


 Milead St Block B  
Faisal Town, Lahore

 0322 - 7561001

 support@medicx4.com

About the company

Our mission is to provide quality health services to all in a comfortable and cozy environment. The clinic primarily offers consultations and is not equipped to take emergencies.



## PART 2 : 醫生功能操作 用姓名查詢

[Home](#) [Patients History](#) [+ Pending Appointments](#) [Today's Appointments](#) [病人基本病歷資料查詢](#) [+ 設定病人量測資料評估數值](#) [病人量測資料查詢](#) [Log Out](#)

姓名  [查詢](#)

	病人姓名	聯絡電話
顯示病人量測資料	ABC	61536516
顯示病人量測資料	DEF	61536516
顯示病人量測資料	XYZ	61536516
顯示病人量測資料	neo	0900654888
顯示病人量測資料	Test0427	0202000000
顯示病人量測資料	陳淳濤	0900666879


姓名: 陳淳濤 [查詢](#)


	病人姓名	聯絡電話
顯示病人量測資料	陳淳濤	0900666879


【醫師量測資料意見評估】  
病人姓名:  
訪查量測資料評估意見

[醫療訊息](#) [最新填單](#)

Healthy Care © 2017


 Milead St Block B  
Faisal Town, Lahore

 0322 - 7561001

 support@medicx4.com

**About the company**

Our mission is to provide quality health services to all in a comfortable and cozy environment. The clinic primarily offers consultations and is not equipped to take emergencies.



## PART 2 : 醫生功能操作 檢視病人歷史量測紀錄資料

備註：  
病人歷史量測紀錄資料  
不會顯示在表中  
請點選【顯示病人量測資料】觀看

姓名  查詢

	病人姓名	聯絡電話
顯示病人量測資料	ABC	61536516
顯示病人量測資料	DEF	61536516
顯示病人量測資料	XYZ	61536516
顯示病人量測資料	neo	0900654888
顯示病人量測資料	Test0427	0202000000
顯示病人量測資料	陳淳潔	0900666879

備註：  
點選【顯示量測時間】  
可觀看填寫時間

姓名  查詢

	病人姓名	聯絡電話
顯示病人量測資料	ABC	61536516
顯示病人量測資料	DEF	61536516
顯示病人量測資料	XYZ	61536516
顯示病人量測資料	neo	0900654888
顯示病人量測資料	Test0427	0202000000
顯示病人量測資料	陳淳潔	0900666879

	No.	量測表單編號	身高	體重	BMI	體溫	心跳脈搏	血氧	血糖	血壓
顯示量測時間	1	6	174	80	26.424	36.5	100	50.78	68.735	95/130
顯示量測時間	2	7	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	3	8	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	4	9	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	5	10	174	90	29.727	36.5	90	100.78	90.735	75/100
顯示量測時間	6	11	180	75	23.148	40	100	50.78	90.735	95/130

	No.	量測表單編號	身高	體重	BMI	體溫	心跳脈搏	血氧	血糖	血壓
顯示量測時間	1	6	174	80	26.424	36.5	100	50.78	68.735	95/130
顯示量測時間	2	7	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	3	8	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	4	9	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	5	10	174	90	29.727	36.5	90	100.78	90.735	75/100
顯示量測時間	6	11	180	75	23.148	40	100	50.78	90.735	95/130

表單填寫時間: 2020/3/5



## PART 2 : 醫生功能操作 給予量測記錄建言



Health Management System

姓名: [Input Field] 全IT

姓名	性別	年齡
陳先生	男	45

No.	量測項目	量測值	單位	標準值	建議
1	體重	80	kg	75	偏高
2	身高	174	cm	174	正常
3	BMI	26.424		24	偏高
4	血壓	120	mmHg	120	正常
5	心跳	78	b/min	78	正常
6	血糖	100	mg/dL	100	正常

【量測結果分析建議】

建議項目: [Input Field]

建議內容: [Input Field]

發送訊息 重新填寫

【量測結果分析建議】

病人姓名: 陳先生

注意紅字警告

發送訊息

重新填寫

## PART 2 : 醫生功能操作 給予量測記錄建言(病人收到的通知)

Home 醫師留言通知 Current Appointment + Take Appointment Treatment History 我的病歷資料 + 填寫量測資料 我的歷史量測資料 Log Out

### 醫師留言通知

There are no new notifications :)

**【醫師叮嚀與囑咐】**  
發言醫師:AJVI  
留言時間:2020/6/2 上午 09:03:54  
留言內容:  
請小心紅字標

**【醫師叮嚀與囑咐】**  
發言醫師:ChenChunYun  
留言時間:2020/6/2 下午 10:02:31  
留言內容:  
注意紅字警告

Health Care © 2017

Milind St Block B,  
Faisal Town, Lahore  
0322 - 7561001  
support@medicx4.com

**About the company**  
Our mission is to provide quality health services to all in a comfortable and cozy environment. The clinic primarily offers consultations and is not equipped to take emergencies.

f

# PART 2 : 醫生功能操作 設定系統量測記錄警訊值

[Home](#) [Patients History](#) [+ Pending Appointments](#) [Today's Appointments](#) [病人基本病歷資料查詢](#) [+ 設定病人量測資料評估數值](#) [病人量測資料查詢](#) [Log Out](#)

[Home](#) [Patients History](#) [+ Pending Appointments](#) [Today's Appointments](#) [病人基本病歷資料查詢](#) [+ 設定病人量測資料評估數值](#) [病人量測資料查詢](#) [Log Out](#)


請填設以下量測數值的警訊範圍

項目	量測數值(最大值)	量測數值(最小值)
體溫	<input type="text" value="請填寫體溫最高值"/>	<input type="text" value="請填寫體溫最低值"/>
脈搏/心跳	<input type="text" value="請填寫脈搏/心跳最高值"/>	<input type="text" value="請填寫脈搏/心跳最低值"/>
血氧	<input type="text" value="請填寫血氧最高值"/>	<input type="text" value="請填寫血氧最低值"/>
血糖	<input type="text" value="請填寫血糖最高值"/>	<input type="text" value="請填寫血糖最低值"/>
血壓(收縮壓)	<input type="text" value="請填寫血壓(收縮壓)最高值"/>	<input type="text" value="請填寫血壓(收縮壓)最低值"/>
血壓(舒張)	<input type="text" value="請填寫血壓(舒張壓)最高值"/>	<input type="text" value="請填寫血壓(舒張壓)最低值"/>


登入設定

If you have any query, please feel free to contact us. ☺


Health Care © 2017.



Milaad St Block B  
Faisal Town, Lahore




0322 - 7661001



support@medicx4.com

About the company

Our mission is to provide quality health services to all in a comfortable and easy environment. The clinic primarily offers consultations and is not equipped to take emergencies.





備註：  
參考資料來源 google 網路資料

## PART 2 : 醫生功能操作 設定系統量測記錄警訊值(範本)

請填設定以下量測數值的警示範圍:

項目	量測數值(最大值)	量測數值(最小值)
體溫	37.9	35.7
脈搏/心跳	100	60
血氧	98	94
血糖	99	60
血壓(收縮壓)	160	140
血壓(舒張)	95	90

## PART 2 : 醫生功能操作

### 設定系統量測記錄警訊值(病人填寫量測記錄提醒)



### 量測記錄表格

請填寫以下量測資料

項目	量測數值	量測時間
身高	請填寫身高(cm)	2020/6/2
體重	請填寫體重(kg)	2020/6/2
BMI BMI參考公式: 體重(公斤) / 身高平方(公尺平方)		2020/6/2
體溫 (標準值: 35.7 ~ 37.9)	請填寫體溫	2020/6/2
脈搏 (標準值: 60 ~ 100)	請填寫脈搏或心跳	2020/6/2
血壓 (標準值: 94 ~ 98)	請填寫血壓	2020/6/2
血糖 (標準值: 60 ~ 99)	請填寫血糖	2020/6/2
血壓收縮壓(高血壓) 收縮(標準值: 140 ~ 160) 舒張(標準值: 90 ~ 95)	請填寫血壓	2020/6/2

送出量測資料

## PART 2 : 醫生功能操作

### 設定系統量測記錄警訊值(病人填寫量測記錄提醒)



量測記錄表單

請填寫以下量測資料

量測日期  
2020/6/2

身高  
請填寫身高(cm)

體重  
請填寫體重(kg)

BMI  
BMI參考公式: 體重(公斤) / 身高平方(公尺平方)

體溫 (標準值: 35.7 ~ 37.9)  
請填寫體溫

脈搏/心跳 (標準值: 60 ~ 100)  
請填寫脈搏或心跳

血氧 (標準值: 94 ~ 98)  
請填寫血氧

血糖 (標準值: 60 ~ 99)  
請填寫血糖

血壓(收縮壓/舒張壓) 收縮 (標準值: 140 ~ 160) 舒張 (標準值: 90 ~ 95)  
請填寫血壓

送出量測資料

## PART 2 : 醫生功能操作

### 設定系統量測記錄警訊值(病人量測記錄警示)



我的歷史量測資料

← → ↻ 不安全 | 182.235.107.91/Patient/PatientMeasurementDataView.aspx

Home 醫師留言通知 Current Appointment Take Appointment Treatment History 我的病歷資料 填寫量測資料 我的歷史量測資料 Log Out

	No.	量測表單編號	身高	體重	BMI	體溫	心跳脈搏	血氧	血糖	血壓
顯示量測時間	1	6	174	80	26.424	36.5	100	50.78	68.735	95/130
顯示量測時間	2	7	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	3	8	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	4	9	174	80	26.424	36.5	120	78.78	35.735	90/110
顯示量測時間	5	10	174	90	29.727	36.5	90	100.78	90.735	75/100
顯示量測時間	6	11	180	75	23.148	40	100	50.78	90.735	95/130

表單填寫時間



## PART 2 : 醫生功能操作 設定系統量測記錄警訊值

---

補充說明：

1. 每個醫生帳號都可以設定系統的警訊值範圍
2. 系統最終採取的警訊值計算方式 -> 所有醫生帳號設定的數值取平均

# 測試站台額外功能 帳號列表

Health Care Login & Registration Form

Health Care Login & Registration Form

Log in to our Website  
Enter Username, Password to log in

Email:  
Password:

OR

Sign up now!  
Fill in the form, create your account

Name:  
Email:  
Address:  
Phone Number:  
Gender:  
Address:  
Password:  
Confirm Password:  
Create Account

or login with:  
Facebook

Search Account Table



System Account Data

Following are the All User Account:

No.	Account	Password	Account Role
1	admin@gmail.com	admin	Admin
2	test1@gmail.com	abc	Doctor
3	test2@gmail.com	abc	Doctor
4	test3@gmail.com	abc	Doctor
5	test4@gmail.com	abc	Doctor
6	test5@gmail.com	abc	Doctor
7	test6@gmail.com	abc	Doctor
8	test7@gmail.com	abc	Doctor
9	test8@gmail.com	abc	Doctor
10	test9@gmail.com	abc	Doctor
11	test10@gmail.com	abc	Doctor
12	test11@gmail.com	abc	Doctor
13	test12@gmail.com	123456	Patient
14	test13@gmail.com	qwe	Patient
15	test14@gmail.com	qwe	Patient
16	test15@gmail.com	123456789	Patient
17	test16@gmail.com	abc	Patient
18	test17@gmail.com	abc	Patient
19	test18@gmail.com	abc	Patient
20	test19@gmail.com	123456789	Patient
21	test20@gmail.com	test123	Patient
22	test21@gmail.com	test123	Patient

以上介紹完畢  
**THANKS**

---

【會後聯絡方式】

E-mail : [neochentwhia@gmail.com](mailto:neochentwhia@gmail.com)